MANUALE D'USO - INSTRUCTION MANUAL

medel[®]

ELITE

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1. INTRODUCTION

The MEDEL ELITE upper arm blood pressure monitor is used to carry out non-invasive measurement and monitoring of the arterial blood pressure values in adults.

This allows you to quickly and easily measure your blood pressure, save the measured values and display the development and average values of the measured values taken. You are also warned of possible existing cardiac arrhythmia. Medel ELITE is a device suitable for home healthcare environment.

ATTENTION!

Please read through this instruction manual carefully, before use. For further questions regarding blood pressure monitoring and measurements, please contact your doctor. Keep Medel ELITE out of the reach of children. Do not use Medel ELITE near anaesthetic, flammable or oxygen mixture, or with noxious oxide.

ATTENTION!

1.1 Important information about self-measurement

- Do not forget: self-measurement means control, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.
- The pulse display is not suitable for checking the frequency of heart pacemakers! In cases of cardiac irregularity (I.H.B.), measurements made with this instrument should only be evaluated after consultation with the doctor.
- Changes to the device are not admitted.
- Do not use the blood pressure monitor on newborns or patients with preeclampsia. We recommend to consult a doctor before using the blood pressure monitor during pregnancy.
- The device does not require any calibration
- Do not dip the device into the water

1.2. Electromagnetic interference

The device contains sensitive electronic components (microcomputer). Therefore, avoid strong electrical or electromagnetic fields close to the device (e.g. mobile phones, microwave ovens...). These can lead to temporary impairment of the accuracy.

2. BLOOD PRESSURE MONITOR COMPONENTS (PIC.1)

- | | |
|---|--------------------------------|
| 1. LCD | 11. Heartbeat symbol |
| 2. O/I button (Start/Stop) | 12. Battery symbol |
| 3. Memory button (MEM) | 13. Systolic blood pressure |
| 4. Average button (T3) | 14. Diastolic blood pressure |
| 5. Air socket | 15. Pulse display |
| 6. Micro USB socket for power | 16. Inflation/Deflation symbol |
| 7. Battery compartment | 17. Date/Time display |
| 8. Universal Cuff (type B applied) for 22-42 cm arm | 18. Stored data |
| 9. Memory symbol | |
| 10. Average value symbol | |

Component list of pressure measuring system:

- PCBA
- Air pipe
- Pump
- Valve
- Cuff

List of accessories:

- Arm cuff
- Storage bag
- 4 x AA Alkaline batteries
- User manual

NOTE:

If your arm cuff is broken or not functional, please use a new cuff. If the new arm cuff does not include an air plug, please continue to use the old air plug.

POWER:

4 AA ALKALINE Batteries (Rechargeable batteries are not suitable for this monitor).

Adaptor (not included) complies with EN 60601.1. Micro USB port for charging*

*  **ATTENTION! THE USB PORT CAN'T BE USED FOR DATA DOWNLOADING!**

Adaptor model not included.

Input: 100-240V~ 50/60Hz 0.4A Max Output: 5V 1000mA

Please only use the authorized adaptor and USB cable to charge power.

Caution: It is not permitted to connect device to any other type of adaptor other than the approved one.

3. PUTTING THE BLOOD PRESSURE MONITOR INTO OPERATION

3.1. Inserting the batteries

Insert the batteries (pic.2)

1. Remove the battery cover from the battery compartment.
2. Insert four AA powerful batteries into the compartment and ensure each battery is in the proper direction.
3. Replace the battery cover.

ATTENTION!

- After the battery warning symbol " " appears, the device is blocked until the batteries have been replaced.
- Please use "AA" long life 1.5 batteries.
- Never leave any low battery into the battery compartment since they may leak and cause damage to the unit.
- If the blood pressure monitor is left unused for long periods, please remove the batteries from the device.
- The typical service life of the new and unused batteries is 450 measurements for the operation time which is 60 s.

3.2 Clock and date adjustment

1. Once you install the battery or turn off the monitor, it will enter Clock mode, and LCD will display date and time.
2. In Clock Mode, press the O/I and MEM buttons simultaneously for 2 seconds and the year value will blink first (Pic.3). Press MEM button to set the year. Keep the

MEM button pressed to increase the numbers.

- Press O/I button to confirm and switch to the month value (press MEM button to set up the month)
- Repeat the procedure described above to set the day, hour and minutes.
- Once you change the batteries, you should readjust the time and date.

3.3 Connecting tube

Insert the cuff connector into the socket on the left side of the meter (Pic.1)

4. HOW TO TAKE PROPER MEASUREMENTS

4.1 Before the measurement:

- Relax and avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors can influence the measurement result.
- Remove any garment that fits closely to your upper arm.
- To take your blood pressure, make sure you are sitting comfortably with your arms and back leaning on something. Do not cross your legs. Place your feet flat on the ground.
- To avoid falsifying the measurement, it is important to remain still during the measurement and not to speak.
- Measure always on the same arm (normally left)
- Make sure to measure the blood pressure always at the same time during the day. It can cause inaccuracy if the measurement is taken in the following circumstances:

- Within 1 hour after dinner or drinking;
- Within 20 minutes after taking a bath;
- When talking or moving your fingers;
- In a very cold environment;

NOTE: The cuff size depends on the circumference of the arm. The allowed range is printed on the cuff. If this is not suitable for your use, please contact your dealer.

⚠ ATTENTION!

- Only use MEDEL original cuffs!
- A loose cuff or a sideways protruding air-pocket causes false measurement values. With repeated measurements, blood accumulates in the respective arm, which can lead to false results.
- Correctly executed blood-pressure measurements should therefore first be repeated after 5 minutes pause or after the arm has been held up in order to allow the accumulated blood to flow away.

4.2. Fitting the cuff

- Pass the end of the cuff (with sewn-in rubber stopper) through the metal stirrup so that a loop is formed. The Velcro closer must be facing outwards. (Ignore this step if the cuff has already been prepared.) (Pic.4)
- Push the cuff over the left upper arm so that the tube points in the direction of the lower arm. (Pic.5)
- Position the cuff on the arm as illustrated in the picture (Pic.6). Make certain that the lower edge of the cuff lies approximately 2 to 3 cm above the elbow and that








the rubber tube leaves the cuff on the inner side of the arm)

- Tighten the free end of the cuff and close the cuff with the closer. (Pic.7)

⚠ ATTENTION!

Place the cuff in order to have its edge 2-3 cm above the elbow and the brachial artery. The tube must be positioned towards the palm center.

- There must be no free space between the arm and the cuff as this would influence the result. Clothing must not restrict the arm. Any piece of clothing which does (e.g. a pullover) must be taken off.
- Secure the cuff with the Velcro closer in such a way that it lies comfortably and is not too tight. Lay the arm on the table (palm upwards) so that the cuff is at the same height as the heart. Make sure that the tube is not kinked. (Pic.8-9).
- Remain seated quietly for two minutes before you begin the measurement.
- The middle of the cuff should be at the level of the right atrium of the heart;
- Before starting measurement, please sit comfortably with legs uncrossed, feet flat on the floor, back and arm supported.

4.3 Measuring procedure

After the cuff has been appropriately positioned, the measurement can begin

1. Press O/I button. All display symbols appear for 1 second. Then 0 will appear on LCD on the DIA column. (pic. 10). The monitor is ready for the measurement.
2. Then the cuff inflates automatically. The symbol "▲" will appear. (pic. 10)
3. Once a pulse is detected, the "♥" symbol blinks. (pic. 11)
4. When appropriate pressure is reached, inflation stops and the cuff starts to deflate.
5. When the measurement is completed, the systolic / diastolic pressure and pulse are shown simultaneously.
6. In case the measurement shows irregular heartbeat (I.H.B.), the "♥" symbol is displayed. In this case the results of the measurement have deviated from a normal situation. Then the measurement needs to be repeated. In the majority of cases, this is not cause for concern. If the symbol "♥" appears regularly (for example several times a week during different daily measurements) it is advisable to discuss it with your medical advisor (pic. 12).

4.4 How to use the memory function:

- In clock mode (Off mode), press MEM button. The number of stored data is shown (Pic.14), sudden the latest result will be displayed with date and time (Pic.15)
- Press MEM button repeatedly to review the results measured.

4.5 T3 System: How to see the average of the last 3 measurements (taken within the last 15 minutes)

1. Press O/I button and take the first measurement.
2. At the end of the first measurement wait at least 5 minutes in a relaxed position. Then repeat the measurement following step 1 again.
3. At the end of the second measurement wait at least 5 minutes in a relaxed position. Then repeat the measurement following step 1 again.
4. At the end of the third measurement press T3 button to perform the average of the last three measurements taken (the AVG symbol will appear) (pic. 16).

NOTE:

The average data using T3 button is shown only if the last 3 measurements have been done within 15 minutes. Otherwise, the data will be shown as 0.

4.6 To delete all values stored in the memory

In Clock mode (Off mode) press MEM button, after the latest data stored will be displayed press again and hold the MEM button for 3 sec. All results will be deleted. Press the MEM or O/I button to turn off the display (pic. 17).

5. ERROR MESSAGE/MALFUNCTION

Follow the below recommended actions and press O/I button again to start another measurement.

Error code	Possible cause(s)	Recommended actions
E1	Deflation is too fast or the cuff doesn't inflate	Re-apply the cuff and try again
E2	The cuff is very tight	
E3	The cuff pressure is excessive	Relax for a moment and try again

E10	The monitor has detected movement during reading	Movement could compromise the reading. Relax for a moment and try again.
E11		
E20	The reading process cannot detect the heart beat	Movement can compromise the measurement. Relax for a moment and try again. Loosen the clothes on the arm and try again
E21	Inaccurate reading	Relax for a moment and try again
E Exx	There is a calibration problem	Repeat the reading. If the problem persists, contact the Service Centre. Refer to the warranty for the contact information and the instructions concerning returns.
Out	Out of measurement range	Relax for a moment. Refasten the cuff and then measure again. If the problem persists, contact your physician.

5.1 Other possible malfunctions:

Problem	Possible cause(s)	Recommended actions
The display shows the symbol 	Batteries are drained.	Replace all batteries with new ones.
The blood pressure values displayed are abnormally high or low.	The cuff is not wrapped correctly or not at the level of the heart.	Wrap the cuff correctly and rise your hand so that the cuff is at the same level of your heart.
	Too much stress is applied on your shoulder or arm.	Relax yourself and make measurement.
	You move your arm or arm muscles during measurement.	Remain still and do not move /contract the muscles during measurement.
The (I.H.B.) symbol "♥" is displayed, but my heart beat rate should be normal.	You move your arm or arm muscles during measurement.	Remain still and do not move /contract the muscles during measurement.
The light on the display is dim or not lit	The batteries are flat or not right placed.	Replace the flat batteries with new ones. Insert the batteries correctly.

6. MAINTENANCE AND REPAIRS

For any repair operations, only refer to a technical assistance centre authorized by the Medel International Srl and request original spare parts. Never open the device. There are no user-serviceable parts within. The device requires no lubrication or maintenance. Failure to comply with what expressed above may compromise the safety of

MEDEL ELITE.

7. CARE, MAINTENANCE AND STORAGE

1. For cleaning only use a soft, lightly moistened cloth.
2. Do not use petroleum spirits, thinners or similar solvents.
3. Do not wash the cuff.
4. Take care to ensure the cuff is not damaged by sharp or pointed objects such as scissors, knives, etc.
5. Take the batteries out if the device is not used for a long period.
6. In case the device becomes malfunction due to the influence of strong electromagnetic fields, power off the device and repeat the procedure again.
7. Disposal of the device (and used batteries) should be carried out according to the national regulations for the disposal of electronic products.
8. Cleaning: Dust environment may affect the performance of the unit. Please use the soft cloth to remove the dirt of the device and cuff before and after use. To obtain the best performance, please follow below instructions:
 - Put in a dry place and avoid the sunshine;
 - Avoid immersing it in the water. Clean it with a dry cloth in case.
 - Avoid shaking and collision.
 - Avoid dusty environment and unstable temperature surrounding.
 - Protect if from:
 - water and moisture
 - extremes temperature, humidity and altitude
 - impact and dropping
 - contamination and dust
 - direct sunlight

8. REFERENCE TO STANDARDS

Device standard: Device corresponds to the requirements of the European standard for non-invasive blood-pressure monitor.

- EN ISO 81060-1:2012 Non-invasive sphygmomanometers - Part 1: Requirements and test methods for non-automated measurement type
- EN 1060-3:1997+A2:2009 Non-invasive blood pressure Part 3: Supplementary requirements for electromechanical blood pressure measuring system
- EN 1060-4: Non invasive sphygmomanometers - Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers.
- IEC 60601-1:2005+A1:2012 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.
- NIBP - Requirements
- IEC 80601-2-30:2009+A1:2013 Medical electrical equipment - Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers
- IEC 60601-1-11:2015 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
- IEC 60601-1-8:2006+A1:2012 Medical electrical equipment - Part 1-8: General

requirements for basic safety and essential performance - Collateral Standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical system.

- IEC 60601-1-2:2007 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard : Electromagnetic compatibility - Requirements and tests

CE 0123 In conformity with the European Medical Device Directive 93/42/EEC

Separate collection for waste electrical and electronic equipment.

Applicable to European Union countries and countries with separate collection systems This symbol on the product or its literature indicates that the product complies with the Electrical and electronic equipment and should not be disposed of as household waste. The user is responsible for the transfer of the equipment at end of life with the appropriate collection facilities, subject to the penalties provided for by current legislation on waste. For more detailed information collection systems available contact your local waste disposal service.

9. PRODUCT SPECIFICATIONS

Display :	LCD large digital backlight display
Technology:	Oscillometric method
Classification :	Internally powered equipment, with type BF applied part
Measurement range:	Rated cuff Pressure : 0 to 300 mmHg
	Measurement Pressure: SYS: 60 mmHg to 230 mmHg DIA: 40 mmHg to 130 mmHg
	Heart rate : 40 to 199 times/minute
Accuracy :	Pressure : +/- 3 mmHg
	Heart rate : +/- 5%
Memory :	Storing 60 memories (Systolic / Diastolic pressure, Heart beat)
Power supply:	4 AA 1.5V (batteries)
Battery life:	For approx. 540 measurements, depending on levels of blood pressure and pump pressure
Accessories:	Universal cuff circum. 22 - 42 cm

Operation conditions:	Temperature : +5°C to +40°C ;
	Relative Humidity: 10 - 85% (non condensing, but not requiring a water vapour partial pressure greater than 50 hPa)
	Pressure: 80 to 105 kPa
Transport/Storage conditions:	Temperature : -20°C to +60°C ;
	Relative Humidity : < 95% (non condensing, at a water vapour pressure up to 50 hPa)
	Pressure: 80 to 105 kPa
Weight of main unit :	390 g (without batteries)
External dimensions:	171 X 131 X 62 mm
Applied part:	Type BF applied part (cuff)
Protection against electric shock:	Internally powered ME equipment /adaptor class II
Protection against harmful ingress of water or particulate matter:	IP22 : The first number 2: protected against solid foreign objects of 12,5 Φ mm and greater. The second number: protected against vertically falling water drops when enclosure tilted up to 15° vertically falling drops.
Mode of operation	Continuous operation
Software version	A7

NOTE:

Specifications are subject to change without prior notice.

Not intended to be sterilized.

Not for use in an OXYGEN RICH ENVIRONMENT.

⚠ ATTENTION!

- To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.
- This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance d away from the equipment. The distance d is calculated by the 800 MHz to 2.5GHz column of Table 6 of IEC 60601-1-2:2007, as appropriate.
- If you have any problems with this device, such as setting up, maintaining or using, please contact with SERVICE PERSONNEL of MEDEL. Don't open or repair the device by yourself. This device must only be serviced, repaired and opened by individuals at authorized sales centers.
- Please report to MEDEL if any unexpected operation or events occur.
- After drop/shock,etc., that may cause changes in the performance, please contact with SERVICE PERSONNEL of MEDEL. Don't open or repair the device by yourself. The use of heat and cold producing devices, such as electric heating blankets,

heating pads or ice packs, may impair the performance of device and increase the risk of injury to the patient.

⚠ ATTENTION!

- Keep the device out of the reach of children/pets to avoid inhalation or swallowing of small parts.
- Be careful to strangulation due to cables and hoses, particularly due to excessive length.
- Do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- Please use ACCESSORIES and detachable parts specified /authorized by MANUFACTURER. Otherwise, it may cause damage to the unit or danger to the user/patients.
- At least 30 min required of equipment to warm from the minimum storage temperature between uses until it is ready for intended use. At least 30 min required for ME equipment to cool from the maximum storage temperature between uses until it is ready for intended use.
- When not in use, store the device with adaptor in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the storage case.
- Do not use solvents to clean this device. Switch off the device and disconnect it from the power supply. The device can be cleaned with a damp cloth. Use lukewarm water and a non-abrasive liquid household cleaner (no abrasive, no alcohol content solution).
- During using, the patient will contact with the cuff. The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010. It will not cause any potential sensitization or irritation reaction.
- Warning: No servicing/maintenance while the ME equipment is in use.
- The patient is an intended operator. The patient can measure data and change battery under normal circumstances and maintain the device and its accessories according to the user manual.
- Adaptor is specified as a part of ME equipment.
- The plug/adaptor plug pins insulates the device from the main supply. Do not position the device in a position where it is difficult to disconnect from the supply mains to safely terminate operation of ME equipment.
- If the pressure of the cuff exceeds 40kPa (300mmHg), the unit will automatically deflate. If the cuff don't deflate when its pressure exceeds 40kPa (300mmHg), detach the cuff from the arm and press the O/I button to stop inflation.
- Before every use, check the device, do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- Warning: No modification of this equipment is allowed.
- Manufacturer will make available on request circuit diagrams, component parts list etc.,.
- The operator shall not touch output of batteries/adaptor and the patient simultaneously.
- The user must check that the equipment does not have visible evidence of damage that may affect patient safety or monitoring capability before use.

- It is intended for adult indoor use only. It is not intended for use with neonatal, pregnant or pre-eclamptic PATIENTS.
- Warning: Be careful regarding the effect of blood flow interference and resulting harmful injury to the patient caused by continuous cuff pressure due to connection tubing kinking; During measurement, please avoid compression or restriction of the connection tubing.
- When using this device, please pay attention to the following situation which may interrupt blood flow and influence blood circulation of the patient, thus cause harmful injury to the patient: connection tubing kinking too frequent and consecutive multiple measurements; the application of the cuff and its pressurization on any arm where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present; inflating the cuff on the side of a mastectomy.
- Warning: Do not apply the cuff over a wound; otherwise it can cause further injury.
- Do not inflate the cuff on the same limb which other monitoring ME equipment is applied around simultaneously, because this could cause temporary loss of function of those simultaneously - used monitoring ME equipment.
- Please check that operation of the device does not result in prolonged impairment of patient blood circulation.
- Do not confuse self-monitoring with self-diagnosis. This unit allows you to monitor your blood pressure. Please start or end medical treatment basing solely on physician's treatment advice. If you are taking medication, consult your physician to determine the most appropriate time for your measurement. Never change a prescribed medication without your physician's consent.
- The device cannot be used with HF surgical equipment at the same time.
- The device is not intended for public use.
- The device is not suitable for continuous monitoring during medical emergencies or operations, as when the cuff is inflated for a long time the patient's blood supply to their arm and fingers will be insufficient, causing anesthesia, distending pain and ecchymosed.
- To verify the calibration of the AUTOMATED SPHYGMOMANOMETER please contact the manufacturer.
- This device may be used only for the purpose described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application.
- The device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in this booklet.
- Do not wash the cuff in a washing machine or dishwasher.
- If Luer lock connectors are used in the construction of tubing, there is a possibility that they might be inadvertently connected to intravascular uid systems, allowing air to be pumped into a blood vessel.
- The patient is an intended operator.
- Expected service life : 2 years
- The service life of the cuff may vary by the frequency of washing, skin condition, and storage state. The typical service life is 10.000 times.

10. INDICATIONS FOR USE

The MEDEL ELITE Blood Pressure Monitor is a digital monitor intended for use in measuring blood pressure and heartbeat rate with arm circumference ranging from 22cm to 42cm (about 8¾"-16½"). It is intended for adult indoor use only.

11. MEASUREMENT PRINCIPLE

This product uses the Oscillometric Measuring method to detect blood pressure. Before every measurement, the unit establishes a "zero pressure" equivalent to the air pressure. Then it starts inflating the arm cuff, meanwhile, the unit detects pressure oscillations generated by beat-to-beat pulsatile, which is used to determine the systolic and diastolic pressure, and also pulse rate. The device also compares the longest and the shortest time intervals of detected pulse waves to mean time interval then calculates standard deviation. The device will display a warning signal with the reading to indicate the detection of irregular heartbeat when the difference of the time intervals is over 25%.

12. SYMBOLS

Refer to the instructions manual

Application part, type BF

For indoor use only

Batch number

CAUTION! These notes must be observed to prevent any damage to the device

Manufacturer's catalogue device number

Serial number

Direct current

Manufacturer

Manufacturer date

Class II Equipment

Keep dry

Temperature limitation for transport and storage

Authorised Representative in the European Community

CE 0123 Complies with the Medical Device Directive 93/42/EEC

Symbol for "ENVIRONMENT PROTECTION - Electrical waste products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice"

13. WARRANTY CONDITIONS

- The device is covered by a 5-year warranty from date of purchase for any defects in the design or in the material used.
- The warranty consists in the replacement and/or repair free of charge of originally defective components.
- The warranty does not cover the accessories supplied with the equipment and those parts subject to normal wear and tear.
- The equipment must be repaired by authorized technical service centres only.
- The equipment must be sent to the authorized technical service centre for repair within 8 days after the defect is noticed.
- The transport costs for the equipment are at user's charge. Any repair out of the warranty conditions are at user's charge.
- The warranty does not cover failures resulting from misuse, abuse of the product or in case the damage is not due to manufacturer's liability (accidental fall, rough transportation, etc...)
- The warranty does not involve any compensation for any kind of damages to persons or things, either direct or indirect, occurred when the equipment is not operating.
- The warranty is valid from the date of purchase certified by the ticket or the invoice which have to be peremptorily attached to the warranty coupon.
- The lack of warranty coupon properly filled in and validated by a certification of purchase voids the warranty.

COUPON TO BE RETURNED IN CASE OF REPAIRS
THE WARRANTY IS VALID ONLY IF ACCOMPANIED BY THE TICKET/INVOICE

Product type: _____

Model: _____

Serial No.: _____

Date of purchase: _____

BUYER'S DATA

Name: _____

Address: _____

Telephone: _____

Defect description: _____

Signature for acceptance of the above mentioned warranty conditions

EMC GUIDANCE

1. This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.
2. * Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
3. Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!
4. * Caution: this machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Table 1) Guidance and MANUFACTURER's declaration
ELECTROMAGNETIC EMISSIONS- for all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration – electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

**Table 2) Guidance and MANUFACTURER's declaration
ELECTROMAGNETIC IMMUNITY – for all ME EQUIPMENT and ME SYSTEMS**

Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below.
The customer or the user of the device should assure that it is used in such an environment

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

**Table 4) Guidance and MANUFACTURER's declaration – electromagnetic IMMUNITY
for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**

Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below.
The customer or the user should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 TEST LEVEL	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. ^a Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Table 6) Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the device.				
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.				
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)			
	150 kHz to 80 MHz $d=1.2 \sqrt{P}$	80 MHz to 800 MHz $d=1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d=2.3 \sqrt{P}$	
0.01	0.12	0.12	0.23	
0.1	0.38	0.38	0.73	
1	1.2	1.2	2.3	
10	3.8	3.8	7.3	
100	12	12	23	
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.				
NOTE 1 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.				
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

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Guangdong Transtek Medical Electronics Co., Ltd.
MDSS - Medical Device Safety Service GmbH